

HELICHECK

Fast Urease Test


The Fastest Method Known

Helicheck
H. Pylori Medium



H. Pylori
NEGATIVE (-)

HP010



Helicheck
H. Pylori Medium



H. Pylori
POZITIVE (+)

• HELICHECK is a fast urease test which is a useful method for detecting *H. Pylori* infections in the endoscopy area. Urease breaks down the urea. The resulting reaction products cause color change in the reagent. Color change is observed between 2 and 30 minutes, depending on the amount of *H. Pylori* the tissue sample. Since the color change varies proportionally depending on the amount of *H. Pylori*, information about the amount of microorganisms in the tissue is obtained.

It is an easy, fast and reliable method for the identification of Helicobacter Pylori

Helicheck is a testing medium that identifies urease activity in gastric biopsy samples

Thanks to its newly developed original formula, it works faster than all other urease tests

It provides optimum conditions in which H. pylori urease is most active and is as sensitive as molecular diagnostic methods